

The background of the entire page is decorated with a network of light blue lines and numerous translucent blue spheres of varying sizes, creating a molecular or network-like aesthetic.

**LEXOGEN**

Enabling complete transcriptome sequencing

**SPLIT**

Fractions for pure RNA sequencing

# RNA Extraction Kit User Guide

Catalog Number:  
008 (SPLIT RNA Extraction Kit)

008UG005V0211

## **FOR RESEARCH USE ONLY. NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE.**

INFORMATION IN THIS DOCUMENT IS SUBJECT TO CHANGE WITHOUT NOTICE.

Lexogen does not assume any responsibility for errors that may appear in this document.

### **PATENTS AND TRADEMARKS**

SPLIT is a trademark of Lexogen. Lexogen is a registered trademark (EU, CH, USA).

Phase Lock Gel is a registered trademark of 5 Prime GmbH, RNAlater is a registered trademark of Ambion, Inc., Austin, Texas, RNasin is a trademark of Promega Corporation, RNaseZap is a registered trademark of Ambion, Inc., Agilent is a registered trademark of Agilent Technologies Inc., Bioanalyzer is a trademark of Agilent Technologies, Inc., DNA-ExitusPlus is a trademark of AppliChem GmbH, Trizol is a registered trademark of Molecular Research Center, Inc. All other brands and names contained in this user guide are the property of their respective owners.

Lexogen does not assume responsibility for violations or patent infringements that may occur with the use of its products.

### **LIABILITY AND LIMITED USE LABEL LICENSE: RESEARCH USE ONLY**

This document is proprietary to Lexogen. The SPLIT RNA Extraction Kits are intended for use in research and development only. They need to be handled by qualified and experienced personnel to ensure safety and proper use. Lexogen does not assume liability for any damage caused by the improper use or the failure to read and explicitly follow this user guide. Furthermore, Lexogen does not assume warranty for merchantability or suitability of the product for a particular purpose.

The purchase of the product does not convey the right to resell, distribute, further sublicense, repackage, or modify the product or any of its components. This document and its contents shall not be used or distributed for any other purpose and/or otherwise communicated, disclosed, or reproduced in any way without the prior written consent of Lexogen.

For information on purchasing additional rights or a license for use other than research, please contact Lexogen.

### **WARRANTY**

Lexogen is committed to providing excellent products. Lexogen warrants that the product performs to the standards described in this user guide until the expiration date. Should this product fail to meet these standards due to any reason other than misuse, improper handling, or storage, Lexogen will replace the product free of charge or issue a credit for the purchase price. Lexogen does not provide any warranty if product components are replaced with substitutes.

Under no circumstances shall the liability of this warranty exceed the purchase price of this product.

### **LITERATURE CITATION**

When describing a procedure for publication using this product, please refer to it as the SPLIT RNA Extraction Kit.

We reserve the right to change, alter, or modify any product without notice to enhance its performance.

---

## **CONTACT INFORMATION**

### **Lexogen GmbH**

Campus Vienna Biocenter 5  
1030 Vienna, Austria  
[www.lexogen.com](http://www.lexogen.com)  
E-mail: [info@lexogen.com](mailto:info@lexogen.com)

### **Support**

E-mail: [support@lexogen.com](mailto:support@lexogen.com)  
Tel. +43 (0) 1 3451212-41  
Fax. +43 (0) 1 3451212-99

# Table of Contents

1. Overview . . . . .	4
2. Kit Components and Storage Conditions . . . . .	6
3. User-supplied Consumables and Equipment. . . . .	7
4. Guidelines . . . . .	8
5. Detailed Protocol . . . . .	10
5.1. Homogenization. . . . .	10
5.1.1. Animal Tissue . . . . .	10
5.1.2. Cell Culture . . . . .	11
5.2. Phenol-Chloroform Extraction . . . . .	12
5.3. Column-based Purification . . . . .	13
5.3.1. Purification of Total RNA . . . . .	13
5.3.2. Purification of Large RNA / Small RNA Fractions . . . . .	15
6. Short Procedure . . . . .	18
7. Appendix A: Tissue Input and Extraction Efficiencies . . . . .	20
8. Appendix B: RNA Quality Control . . . . .	20
9. Appendix C: RNA Storage . . . . .	22
10. Appendix D: Revision History . . . . .	23

# 1. Overview

The SPLIT RNA Extraction Kit enables a fast and highly efficient extraction of RNA that is free of genomic DNA contamination. The RNA can be recovered as total RNA or split into a large and a small RNA fraction, facilitating the analysis of e.g., mRNA and miRNA from the same sample. Thus the RNA obtained is ideal for seamlessly preparing libraries for Next Generation Sequencing of total RNA or its large and small fractions. Importantly, the SPLIT protocol does not require DNase treatment for the removal of genomic DNA and thereby avoids the high risk of concomitant RNA degradation.

First, the sample is homogenized in an isolation buffer that is highly chaotropic to facilitate effortless and complete solubilization.

Acidic buffer and acidic phenol are added to create a monophasic solution, a step that is essential for the efficient separation of genomic DNA into the organic phase. Chloroform is added, and phases are cleanly separated using phase lock gel tubes. The use of these tubes mitigates the risk of contaminating the upper aqueous phase that contains RNA with the lower phenol phase that contains DNA and protein.

The RNA is further purified on a silica column to eliminate trace amounts of phenol and to optionally fractionate the RNA. By adding 1.75 x volume of isopropanol to the aqueous phase the entire total RNA will precipitate onto the silica carrier. When using only 0.33 x volume isopropanol large RNA with a lower limit of about 150 nt will bind whereas the small RNA will be in the flow-through. By adding 1x volume of isopropanol to this flow-through also the small RNA can be recovered on another silica column.

With the SPLIT RNA Extraction Kit, either 48 samples can be extracted for their total RNA or large RNA fraction only (single-fraction extractions), or small and large RNA fractions can be isolated from 24 samples (dual-fraction extractions).

Protocols are given for RNA extraction from human cell culture and animal tissue, and the extraction protocol can be easily adapted to a variety of RNA sources (Appendix A, p.20). Please note that acidic phenol, chloroform, ethanol, and isopropanol have to be supplied by the user.

The RNA obtained is of highest purity and ideally suited to prepare libraries for RNA sequencing. It is also superior for other high demanding applications such as full-length reverse transcription or sample preparation for microarray analysis (see Appendix B, p.20 for details on RNA quality).




Figure 1. Schematic overview of the SPLIT workflow. The cut-off between large and small RNA fractions is at 150 nt.

## 2. Kit Components and Storage Conditions




Figure 2. Location of kit contents.

Kit Component	Label	Requirement for		Storage
		48 single-fraction purifications	24 dual-fraction purifications	
Isolation Buffer	IB	19.2 ml	9.6 ml	+4 °C
Acidic Buffer	AB	7.2 ml	3.6 ml	+4 °C
Wash Buffer	WB	3 x 17.6 ml <sup>1</sup>	3 x 17.6 ml <sup>1</sup>	+4 °C
Elution Buffer	EB	2.4 ml <sup>2</sup>	2.4 ml <sup>2</sup>	+4 °C
Storage Buffer	SB	2.4 ml <sup>2</sup>	2.4 ml <sup>2</sup>	+4 °C
Phase lock gel tubes	Phase lock gel tubes	48	24	+4 °C
Purification columns	Purification columns	48	48	+4 °C

<sup>1</sup> Including ethanol added by the user.

<sup>2</sup> For each RNA fraction, either EB or SB is required.

Buffers provided in the SPLIT kit include tolerance volumes.

Upon receiving the SPLIT kit, store it in a +2 to +8 °C refrigerator.

**Caution:** Phase lock gel tubes must not be frozen.

**IB** is to be used at +4 °C. All other components (especially phase lock gel tubes) should get to room temperature before use.

Check the contents of **IB**, **AB**, **WB**, and **SB** which may precipitate during shipping and storage. If a white precipitate is visible, incubate at 37 °C until buffer components dissolve completely.

Cat. No. 008.48 (48 extractions): Add 20 ml absolute ethanol to each of the 3 bottles with Wash Buffer (**WB**) concentrate and shake to combine. This will yield 3 x 25 ml (75 ml) Wash Buffer (**WB**).

## 3. User-supplied Consumables and Equipment

Check to ensure that you have all of the necessary materials and equipment before beginning with the RNA extraction. All reagents, equipment and labware must be free of nucleases and nucleic acid contamination.

### Reagents

Reagent	Requirement for			Comment
	24 dual fraction purifications	48 large RNA only purifications	48 total RNA purifications	
<b>Phenol solution pH 4.3</b>	9.6 ml	19.2 ml	19.2 ml	e.g. Sigma-Aldrich P4682-100ML
<b>Chloroform</b>	4.8 ml	9.6 ml	9.6 ml	
<b>Isopropanol</b>	24.0 ml	9.6 ml	50.4 ml	2-Propanol
<b>Ethanol abs.</b>	3 x 20 ml	3 x 20 ml	3 x 20 ml	Added to <b>WB</b>

### Equipment

- Fume hood for organic solvent handling.
- Benchtop centrifuge (12,000 x g, rotor compatible with 1.5 ml and 2.0 ml micro-tubes).
- Calibrated single-channel pipettes for handling 10 µl to 1000 µl volumes.
- Vortex mixer.
- UV-spectrophotometer to quantify RNA.

### Labware

- Suitable pipette tips (pipette tips with aerosol barriers recommended).
- 1.5 ml and 2.0 ml micro-tubes with cap, low binding, certified ribonuclease-free.
- Benchtop cooler or ice pellets in ice box (for short-term storage of RNA).

### Optional Equipment & Solutions

- 0.1 ml tissue grinder (hand-held homogenizer).
- Automated microfluidic electrophoresis station (Agilent Technologies 2100 Bioanalyzer).
- Agarose gels, dyes, and electrophoresis rig (for RNA quality control).
- DNA-ExitusPlus (AppliChem GmbH).
- RNaseZap.
- RNase inhibitor.

The complete set of materials, reagents, and labware for quality control is not listed. Consult 5.1. (p.10) for homogenization protocols and Appendix B (p.20) for more information on RNA quality.

## 4. Guidelines

### Safety Information

- This kit is to be used with a phenol solution, which is toxic and corrosive and with chloroform, both of which may be health hazards if not handled properly. Phenol should not come in contact with skin, eyes, or the respiratory tract and may cause chemical burns to the exposed area. When working with the phenol solution and with chloroform, always work in a fume hood.
- The Isolation Buffer (**IB**) contains guanidine isothiocyanate, an irritant, which upon protocol completion is also present in flow-through and wash fractions. This chemical is harmful when in contact with the skin, inhaled, or ingested. Do not add bleach or acidic solutions directly to solutions or sample preparation waste that contains guanidine isothiocyanate, as reactive compounds and toxic gases are formed.
- Solutions containing isopropanol or ethanol are considered flammable. Use appropriate precautions when using these chemicals.
- For your protection, always wear a laboratory coat, gloves, and safety glasses when handling chemicals.
- Dispose of the buffers and chemicals in appropriate waste containers.
- Consult the appropriate Material Safety Data Sheets (MSDS), available at [www.lexogen.com](http://www.lexogen.com), and contact your Environmental Health and Safety department for proper work and disposal guidelines.

### General

- Unless explicitly mentioned, all steps should be carried out at a room temperature (RT) between 20 °C and 25 °C. Centrifugation should be performed at 18 °C to increase reproducibility. If a refrigerated centrifuge is not available, centrifugation can be carried out at RT. Results may be negatively impacted if the protocol is performed at temperatures outside of this range.
- Ensure that adequate volumes of all reagents and the necessary equipment are available before beginning the protocol.
- Perform all pipetting steps with calibrated pipettes and always use fresh tips. Pipette carefully to avoid foaming as the Isolation Buffer (**IB**) and the Storage Buffer (**SB**) contain detergents.
- Before you start, check all solutions for the formation of precipitate and if necessary, incubate at 37 °C until buffer components dissolve completely.
- The phenol solution pH 4.3 and the Isolation Buffer (**IB**) should be used at their storage temperature of +4 °C. All other components, especially the phase lock gel tubes, should be at room temperature.



## RNA Handling

- RNases are ubiquitous, and special care should be taken throughout the procedure to avoid RNase contamination.
- Use commercial ribonuclease inhibitors (i.e., RNasin, Promega Corp.) to maintain RNA integrity when storing samples.
- Use a sterile and RNase-free workstation or laminar-flow hood if available. Please note that RNases may still be present on sterile surfaces, and that autoclaving does not completely eliminate RNase contamination.
- Before starting an RNA extraction, clean your work space, pipettes, and other equipment with RNase removal spray (such as RNaseZap, Ambion Inc.) as per the manufacturer's instructions.
- Protect all reagents and your RNA samples from RNases on your skin by wearing a clean lab coat and fresh gloves. Change gloves after making contact with equipment or surfaces outside of the RNase-free zone.
- Avoid speaking above opened tubes. Keep reagents closed when not in use to avoid airborne RNase contamination.

# 5. Detailed Protocol

## 5.1. Homogenization

### 5.1.1. Animal Tissue

#### Preparation

Tissue	Weigh and Reduce Tissue	Homogenization
<b>Animal tissue</b> – freshly harvested or frozen or thawed at +4 °C if stored in RNAlater	<b>Tweezers</b> – sterile <b>Scalpel</b> – sterile <b>Gauze pad</b> – sterile	<b>Isolation Buffer (IB)</b> - at +4 °C or on ice
<b>Fume hood</b> or laminar-flow cabinet	<b>Precision balance</b>	<b>Tissue grinder</b> – 0.1 ml, glass

#### Homogenization

Tissue is homogenized in a highly chaotropic solution.

This protocol is specific for hand-held tissue grinders (glass homogenizers with pestle) but can be easily adapted for other homogenization protocols. Optimally, the tissue should be stored at -20 °C in RNAlater (Ambion Inc.). Tissue frozen without preservation (RNAlater) must not be thawed before homogenization to maintain RNA integrity. To prevent cross-contamination, it is best to work in a fume hood or a laminar-flow cabinet that can be UV-irradiated.

- 1 Add 400 µl cold (+4 °C) Isolation Buffer (**IB**) into a glass tissue grinder.
- 2 Use sterile tweezers to transfer a tissue piece onto a fresh, sterile gauze pad. If RNAlater was used for conservation, dry the tissue by tapping onto the gauze pad.
- 3 Determine the weight of the tissue on a precision balance. The protocol is efficient for extraction of up to 100 µg of total RNA. See Appendix A (p.20) for details on input and extraction efficiency.
- 4 Optional: Hard to homogenize tissue such as tendons or cartilage can be reduced using a scalpel to facilitate solubilization in the next steps.
- 5 Using tweezers, transfer the tissue pieces quantitatively into the Isolation Buffer (**IB**) in the tissue grinder.
- 6 Homogenize the tissue by carefully moving the pestle up and down. Simultaneous rotation helps to dissolve also larger pieces. Do not pull out the pestle completely to avoid foaming. The tissue is usually homogenized within 2 - 3 minutes; avoid extended homogenization and warming up of Isolation Buffer (**IB**).

- 
- 7 Continue immediately with the phenol-chloroform extraction in 5.2.
- 

After use, clean the tissue grinder thoroughly with a detergent such as DNA-ExitusPlus (Appli-Chem GmbH), then with ultra-filtered water and finally with 75% ethanol.

## 5.1.2. Cell Culture

### Preparation

Cells	Solubilization
Cells – freshly harvested or frozen or thawed at +4 °C if stored in RNeasy lysis buffer	Isolation Buffer (IB) – at +4 °C or on ice
Fume hood or laminar-flow cabinet	

### Solubilization

Cells are solubilized in a highly chaotropic solution.

If cells are not harvested freshly, they can be pelleted, washed and stored in RNeasy lysis buffer (Ambion, Inc.) at -20 °C. To prevent cross-contamination, it is best to work in a fume hood or a laminar-flow cabinet that can be UV-irradiated.

- 1 Harvest, pellet and wash the cells. If they have been stored in RNA later, they should be pelleted at +4 °C and RNeasy lysis buffer removed with a pipette. The protocol is suitable for extraction of e.g.,  $10^6$  cells of a human suspension cell culture.
  - 2 Add 400 µl cold (+4 °C) Isolation Buffer (IB) to the cells.
  - 3 Lyse the cells by carefully pipetting up and down. The cells are usually lysed within 1 - 2 minutes.
  - 4 Continue with step 8, phenol-chloroform extraction in 5.2.
-

## 5.2. Phenol-Chloroform Extraction

### Preparation

	For each sample	Temperature
Phenol solution pH 4.3 <sup>1</sup>	400 µl	+4 °C
Acidic buffer (AB)	150 µl	RT
Chloroform <sup>1</sup>	200 µl	RT
Phase lock gel tube	1	RT
2 ml micro-tube	1	RT
Centrifuge		18 °C
Fume hood		
Vortex mixer		

<sup>1</sup> **Caution:** When working with phenol or chloroform always use a fume hood and discard waste according to applicable Health and Safety regulations.

### Phenol-Chloroform Extraction

Utilizing a highly specific phenol-chloroform extraction, RNA is partitioned into the upper, aqueous phase whereas DNA and proteins are partitioned into the lower, organic phase.

- 8 For each sample, centrifuge 1 phase lock gel tube for 1 minute at 12,000 x g at 18 °C. This collects the gel on the bottom of the tube.
- 9 Transfer the homogenized tissue (cells) in Isolation Buffer (**IB**) into a phase lock gel tube.
- 10 Add 400 µl phenol solution pH 4.3 and mix by inverting the tube 5 times.
- 11 Add 150 µl Acidic Buffer (**AB**) and mix by pipetting.
- 12 Add 200 µl of chloroform.  
  
Mix vigorously by 3 cycles of 5 seconds vortexing and 1 second pausing.  
**Attention:** Vigorous vortexing is essential to disperse the chloroform efficiently and effectively separate all the phenol that will contain the gDNA and protein into the organic phase. Do not be afraid of shearing the gDNA. Even if this happens all DNA will separate into the lower organic phase irrespective of its size.
- 14 Incubate for 2 minutes at room temperature.
- 15 Centrifuge for 2 minutes at 12,000 x g at 18 °C.
- 16 Transfer the upper phase to a new 2 ml micro-tube by decanting.  
  
For the purification of **total RNA**, proceed with **step 18 in 5.3.1**. For the purification of the **large RNA fraction** and – optionally – the **small RNA fraction**, proceed with **step 18 in 5.3.2**.

## 5.3. Column-based Purification

### 5.3.1. Purification of Total RNA

#### Preparation

	For each sample	Temperature
Isopropanol	~1050 µl	RT
Wash Buffer (WB) <sup>1</sup>	1100 µl	RT
Elution Buffer (EB) or Storage Buffer (SB) <sup>2</sup>	50 µl	RT
Purification column	1	RT
Collection tube	1	RT
1.5 ml micro-tube	1	RT
Centrifuge		18 °C
Vortex mixer		

<sup>1</sup> **Caution:** Discard waste containing guanidine isothiocyanate, phenol and chloroform according to applicable Health and Safety regulations.

<sup>2</sup> see Appendix C (p.22) whether **EB** or **SB** should be used for elution.

#### Column Loading

The total RNA is precipitated onto a silica column by addition of 1.75 x volume of isopropanol.

- 18 Determine the volume of the transferred upper phase in the 2 ml micro-tube and add isopropanol at 1.75 x of this volume (e.g. 1050 µl isopropanol to 600 µl sample). Mix by vortexing for 10 seconds.
- 19 Place a purification column in a collection tube.
- 20 Apply a maximum of 800 µl of the mixture from step 18 (aqueous phase with isopropanol) to the column.
- 21 Centrifuge for 20 seconds at 12,000 x g at 18 °C and discard the content of the collection tube.
- 22 Repeat steps 20 - 21 until the mixture is loaded completely.

## Column Washing

The RNA is further purified by washing on the column.

- 23 Apply 500 µl of Wash Buffer (**WB**) to the column and centrifuge for 20 seconds at 12,000 x g at 18 °C. Empty the collection tube.
- 24 Apply 600 µl of Wash Buffer (**WB**) to the column and centrifuge for 20 seconds at 12,000 x g at 18 °C.
- 25 Empty the collection tube with a pipette and centrifuge for 1 minute at 12,000 x g at 18 °C.
- 26 Discard the collection tube and place the purification column in new 1.5 ml micro-tube. Make sure that no ethanol traces are carried to the new tube.

## Elution of Total RNA

The RNA is eluted into an elution or storage buffer.

- 27 Apply either 10 - 50 µl Elution Buffer (**EB**) or 10 - 50 µl Storage Buffer (**SB**) to the column. Incubate for 1 minute at room temperature.
- 28 Centrifuge for 1 minute at 12,000 x g at 18 °C.
- 29 Optional: For a second elution, place the purification column in a new 1.5 ml micro-tube and repeat steps 27 - 28.
- 30 At this point the total RNA is purified and ready for quality control (Appendix B, p.20) and downstream applications.
- 31 Optional: Add RNase inhibitor (not included). See Appendix C (p.22) for RNA storage. Note that the RNase inhibitor might absorb at 230 nm, therefore use buffer with RNase inhibitor added as blank in OD measurements.

## 5.3.2. Purification of Large RNA / Small RNA Fractions

### Preparation

	Large RNA fraction	Large RNA and small RNA fraction	Temperature
<b>Isopropanol</b>	~200 µl	~800 µl	RT
<b>Wash Buffer (WB)<sup>1</sup></b>	1100 µl	2200 µl	RT
<b>Elution Buffer (EB) or Storage Buffer (SB)<sup>2</sup></b>	(2x) 50 µl	(2x) 100 µl	RT
<b>Purification column(s)</b>	1	2	RT
<b>Collection tube(s)</b>	1	2	RT
<b>2 ml micro-tube</b>	-	1	RT
<b>1.5 ml micro-tube(s)</b>	1	2	RT
<b>Centrifuge</b>			18 °C
<b>Vortex Mixer</b>			

<sup>1</sup> **Caution:** Discard waste containing guanidine isothiocyanate, phenol and chloroform according to applicable Health and Safety regulations.

<sup>2</sup> see Appendix C (p.22) whether **EB** or **SB** should be used for elution.

### Column Loading

The large RNA fraction is precipitated onto a silica column by the addition of 0.33 x volume of isopropanol. The small RNA fraction will be in the flow-through.

- 18 Determine the volume of the aqueous phase, which may vary, depending on the tissue volume and volume transfer efficiency during homogenization and extraction.  
**Attention:** For best reproducibility of the size cut-off it is essential to quantify the volume of the aqueous phase exactly.

- 19 Add isopropanol at 0.33 x of this volume (e.g. 200 µl isopropanol to 600 µl sample). Mix by vortexing for 10 seconds.

- 20 Place a purification column in a collection tube.

- 21 Apply a maximum of 800 µl of the mixture (aqueous phase / isopropanol) to the column.

- 22 Centrifuge for 20 seconds at 12,000 x g at 18 °C.

- 23 **Optional: If you want to also isolate the small RNA fraction pipette the flow-through into a 2 ml micro-tube. Else discard the flow-through.**

- 24 Repeat steps 21 - 23 until the mixture is loaded completely.

## Splitting into Fractions

The small RNA fraction is recovered by precipitation onto a new purification column with the addition of 1 x volume of isopropanol.

- 25** **LARGE RNA ONLY:** If you want to only isolate the large RNA fraction, then you can discard the content of the collection tube. Place the purification column back into the collection tube and proceed with step 28.
- 

- 26** **SMALL RNA PURIFICATION (optional).** If you want to also isolate the small RNA fraction; then measure the volume of the flow-through in the 2 ml micro-tube and add the same volume isopropanol (e.g. 800 µl isopropanol to 800 µl flow-through). Mix by vortexing for 10 seconds.
- 

- 27** **SMALL RNA PURIFICATION (optional).** For loading of the small RNA fraction, place a new purification column in a collection tube. Apply a maximum of 800 µl of the flow-through / isopropanol mixture to the column. Centrifuge for 20 seconds at 12,000 x g at 18 °C and discard the content of the collection tube. Repeat until loading of the mixture is complete. Washing and elution can now be performed for the large RNA fraction and the small RNA fraction in parallel.
- 

## Column Washing

The RNA fractions are further purified by washing on the column.

This protocol applies to both, the large RNA fraction on the original purification column and – optionally – the small RNA fraction on the second purification column.

- 28** Apply 500 µl of Wash Buffer (**WB**) to the purification column and centrifuge for 20 seconds at 12,000 x g at 18 °C. Empty the collection tube.
- 
- 29** Apply 600 µl of Wash Buffer (**WB**) to the column and centrifuge for 20 seconds at 12,000 x g at 18 °C.
- 
- 30** Empty the collection tube with a pipette and centrifuge for 1 minute at 12,000 x g at 18 °C.
- 
- 31** Discard the collection tube and place the purification column in new 1.5 ml micro-tube. Make sure that no ethanol traces are carried to the new tube.
-



## Elution of RNA Fractions

The RNA fractions are eluted into an elution or storage buffer.

- 32 Apply either 10 - 50  $\mu$ l Elution Buffer (**EB**) or 10 - 50  $\mu$ l Storage Buffer (**SB**) to each column. Incubate for 1 minute at room temperature.

---

- 33 Centrifuge for 1 minute at 12,000 x g at 18 °C.

---

- 34 Optional: For a second elution, place the purification column in a new 1.5 ml micro-tube and repeat steps 32 - 33.

---

- 35 At this point the large RNA fraction – and optionally, also the small RNA fraction – is purified and ready for quality control (Appendix B, p.20) and downstream applications.

---

- 36 Optional: Add RNase inhibitor (not included). See Appendix C (p.22) for RNA storage. Note that the RNase inhibitor might absorb at 230 nm, therefore use buffer with RNase inhibitor added as blank in OD measurements.

---

## 6. Short Procedure

### 6.1. Extraction of Total RNA

All centrifugation steps are at 12,000 x g and 18 °C.

#### 20 min Homogenization and Phenol-Chloroform Extraction

<input type="checkbox"/> Homogenize tissue / cells in 400 µl <b>IB</b> .	Homogenization
<input type="checkbox"/> Centrifuge 1 phase lock gel tube for 1 min.	
<input type="checkbox"/> Transfer homogenate into a phase lock gel tube.	
<input type="checkbox"/> Add 400 µl phenol solution pH 4.3, mix by inverting the tube 5 times.	
<input type="checkbox"/> Add 150 µl <b>AB</b> , mix by pipetting.	Phenol-Chloroform Extraction
<input type="checkbox"/> Add 200 µl chloroform and mix vigorously (3x 5 sec vortexing and 1 sec pausing).	
<input type="checkbox"/> Incubate for 2 min at RT.	
<input type="checkbox"/> Centrifuge for 2 min.	
<input type="checkbox"/> Decant the upper phase into a 2 ml micro-tube.	

#### 10 min

#### Purification of Total RNA

<input type="checkbox"/> Measure volume of transferred upper phase.	Column Loading
<input type="checkbox"/> Add 1.75x vol. isopropanol to the upper phase.	
<input type="checkbox"/> Mix by vortexing for 10 sec.	
<input type="checkbox"/> Load max. 800 µl onto purification column in collection tube.	
<input type="checkbox"/> Centrifuge for 20 sec and discard flow-through. Repeat until mixture is loaded completely.	
<input type="checkbox"/> Apply 500 µl <b>WB</b> and centrifuge for 20 sec. Empty collection tube.	Column Washing
<input type="checkbox"/> Apply 600 µl <b>WB</b> and centrifuge for 20 sec.	
<input type="checkbox"/> Empty collection tube (with pipette) and centrifuge for 1 min.	
<input type="checkbox"/> Place purification column in a 1.5 ml micro-tube.	Elution
<input type="checkbox"/> Apply 10 - 50 µl of either <b>EB</b> or <b>SB</b> , incubate for 1 min at RT.	
<input type="checkbox"/> Centrifuge for 1 min.	
<input type="checkbox"/> Optional: Second elution in new 1.5 ml micro-tube.	
<input type="checkbox"/> Optional: Add RNase inhibitor (not included).	

# 6.2. Extraction of Large RNA Fraction or Large & Small RNA Fractions

All centrifugation steps are at 12,000 x g and 18 °C.

## 20 min Homogenization and Phenol-Chloroform Extraction

<input type="checkbox"/>	Homogenize tissue / cells in 400 µl <b>IB</b> .	Homogenization
<input type="checkbox"/>	Centrifuge 1 phase lock gel tube for 1 min.	
<input type="checkbox"/>	Transfer homogenate into a phase lock gel tube.	
<input type="checkbox"/>	Add 400 µl phenol solution pH 4.3, mix by inverting the tube 5 times.	
<input type="checkbox"/>	Add 150 µl <b>AB</b> , mix by pipetting.	Phenol-Chloroform Extraction
<input type="checkbox"/>	Add 200 µl chloroform and mix vigorously (3x 5 sec vortexing and 1 sec pausing).	
<input type="checkbox"/>	Incubate for 2 min at RT.	
<input type="checkbox"/>	Centrifuge for 2 min.	
<input type="checkbox"/>	Decant the upper phase into a 2 ml micro-tube.	

## 15 min Purification of Large / Small RNA Fraction(s)

<input type="checkbox"/>	Measure volume of transferred upper phase.	
<input type="checkbox"/>	Add 0.33x vol. isopropanol to the upper phase.	
<input type="checkbox"/>	Mix by vortexing for 10 sec.	Column Loading
<input type="checkbox"/>	Load max. 800 µl onto purification column in collection tube.	
<input type="checkbox"/>	Centrifuge for 20 sec. Transfer flow-through into a 2 ml micro-tube. Repeat until loading is complete.	
LARGE RNA FRACTION ONLY: Discard 2 ml micro-tube.		
ADDITIONAL SMALL RNA PURIFICATION:		
<input type="checkbox"/>	<input type="checkbox"/> Measure flow-through volume in 2 ml micro-tube and add 1x vol. isopropanol. Mix by vortexing for 10 seconds.	Fraction Split
<input type="checkbox"/>	<input type="checkbox"/> Load max. 800 µl onto (new) purification column in a collection tube. Centrifuge for 20 sec. Discard flow-through and repeat until loading is complete.	
<input type="checkbox"/>	Apply 500 µl <b>WB</b> to each column and centrifuge for 20 sec. Empty collection tube(s).	Column Washing
<input type="checkbox"/>	Apply 600 µl <b>WB</b> to each column and centrifuge for 20 sec.	
<input type="checkbox"/>	Empty collection tube(s) (with pipette) and centrifuge for 1 min.	
<input type="checkbox"/>	Place purification column(s) in 1.5 ml micro-tube(s).	
<input type="checkbox"/>	Apply 10 - 50 µl of either <b>EB</b> or <b>SB</b> , incubate for 1 min at RT.	
<input type="checkbox"/>	Centrifuge for 1 min.	Elution
<input type="checkbox"/>	Optional: Second elution in new 1.5 ml micro-tube.	
<input type="checkbox"/>	Optional: Add RNase inhibitor (not included).	

## 7. Appendix A: Tissue Input and Extraction Efficiencies

If immediate RNA extraction is not possible, tissue samples can be either flash-frozen with liquid nitrogen and stored at -80 °C or preserved in RNAlater (Ambion, Inc.) and stored at -20 °C or -80 °C. Tissue / cells without RNAlater preservation must only be thawed during the homogenization step in cold Isolation Buffer (+4 °C) to keep RNases inactive.

RNA extraction efficiency for mouse liver is typically 4.0 - 4.5 µg total RNA / mg tissue (3.0 - 3.5 µg large RNA and 0.6 µg small RNA / mg tissue). A second elution can increase the overall yield by 5 - 20 %. The maximum binding capacity of the purification column is 100 µg RNA, which should not be exceeded for optimal results. For mouse liver tissue this translates into an upper limit of 20 - 25 mg input per extraction. Other tissues have different RNA content, and the input might have to be adjusted accordingly.

## 8. Appendix B: RNA Quality Control

### RNA Integrity

The integrity of an RNA sample can be assessed with a variety of methods (see Figure 3). We recommend the use of a microfluidics assay such as the RNA6000 series for the 2100 Bioanalyzer (Agilent Technologies Inc.), although RNA quality can also be assessed with denaturing agarose gel electrophoresis if such a device is not available. Most microfluidics platforms will carry out an automated peak analysis and generate a quality score (RIN or RQN) in addition to the 28S/18S rRNA ratio. The quality of RNA extracted with the SPLIT RNA Extraction Kit almost exclusively depends on the extraction source: a RIN of 10 and a 28S/18S rRNA ratio of 2.7 can be obtained from human cell culture homogenized according to 5.1.2. Extractions from tissue samples usually result in RNA with a RIN of 8.0 - 9.5.


### Potential Contaminants

RNA samples should be free of salts, metal ions, and organic solvents, which can be carried over from the RNA extraction. Several sources of contamination can be detected with a UV-Vis spectrophotometer. An acceptably pure RNA sample should have an A260/A280 ratio between 1.8 and 2.1. The A260/A230 ratio should also be greater than 1.8. Several common contaminants including proteins, chaotropic salts and phenol absorb strongly between 220 and 230 nm and can often be identified as peaks in this region. Contamination with any of these generates a lower A260/230 ratio. Phenol also has an absorption maximum between 250 and 280 nm, which overlaps that of nucleic acid, so high 230 nm absorbance combined with a biphasic or broad peak between 250 and 280 nm may indicate contamination with phenol rather than chaotropic salts.

# Genomic DNA Contamination

The SPLIT RNA Extraction Kit was designed for minimizing the genomic DNA (gDNA) content in the RNA sample. gDNA is indistinguishable from RNA on a spectrophotometer, and many of the dyes used in RNA microfluidics assays stain single-stranded nucleic acids much more intensely than double-stranded. Hence, low to moderate amounts of gDNA may not be readily visible with an RNA-specific microfluidics assay. We highly recommend examining all RNA samples on a denaturing agarose gel or using a fluorometric assay with DNA- and RNA-specific dyes to check samples for DNA contamination. On an agarose gel, gDNA can appear as either a dark mass which remains in the slot if relatively intact (see Figure 3B) or as a high molecular weight smear if it has been sheared during extraction.

## Typical Results



E	Tissue	RNA yield	Efficiency*	RIN	A260/A280	A260/A230
Total RNA	19.8 mg	88.0 µg	4.43 µg/mg	8.2	1.98	2.22
Large RNA	19.8 mg	64.7 µg	3.26 µg/mg	8.3	2.00	2.28
Small RNA	19.8 mg	12.3 µg	0.62 µg/mg	N/A	2.00	2.06

**Figure 3. Analysis of SPLIT kit extracted RNA.** (A) Gel-like representation of Agilent Bioanalyzer traces. RNA from mouse liver stored in RNAlater was extracted either as total RNA (lane 1) or as large RNA and small RNA fractions (lanes 2 and 3). In the split sample RNAs shorter than 150 nt are confined to the small RNA fraction. A control sample was extracted following a TRIzol protocol (lane 4). This RNA sample contains a significant amount of genomic DNA (gDNA) that is not detected by the Bioanalyzer but becomes visible on a denaturing agarose gel as a slot-retained band. (B) RNA obtained with the SPLIT kit is free from detectable gDNA. (C) A miRNA marker was spiked into mouse liver homogenate, which was then extracted using the SPLIT kit. Analysis on a 15 % denaturing polyacrylamide gel demonstrates that small RNA down to at least 17 nt is efficiently recovered in the total RNA sample and in the small RNA fraction. The theoretical maximum spike-in RNA recovery amount was loaded in lane 1. (D) Bioanalyzer evaluation of miRNA-spiked samples on a small RNA chip (10 - 200 nt, linear scale) and on an RNA 6000 pico chip (200 - 500 nt, log scale). The traces from the two chips are shown alongside for illustrative purposes, the Y-axes do not correspond quantitatively. (E) The table shows key parameters of SPLIT RNA extractions from mouse liver that was stored in RNAlater. Efficiency is given in µg RNA per mg tissue.

## 9. Appendix C: RNA Storage

After extraction, RNA can be stored in Elution Buffer (**EB**, 10 mM Tris-HCl pH 7.0) at -20 °C or -80 °C. This minimal buffer stabilizes the pH without any other components that might interfere with downstream applications. When eluting in **EB** we highly recommend the addition of RNase inhibitor to block any accidentally introduced RNases.

The Storage Buffer (**SB**, 10 mM Tris-HCl pH 7.0, 10 mM DTT and 0.1 mM EDTA) supplied with this kit can be used for intermediate storage of the RNA at -20 °C or -80 °C. DTT (antioxidant) and EDTA (chelating agent) both minimize the threat of RNA degradation, especially at non-freezing conditions. For long-term storage, we recommend to keep aliquots of the RNA as NaAc/ethanol precipitate at -80 °C to avoid accidental RNase contamination as well as RNA degradation due to freeze/thaw cycles.

We suggest to check the RNA quality after extended periods of storage for changes in integrity and quantity e.g., on a microfluidics system.

# 10. Appendix D: Revision History

Revision date	Change	Page(s)
July 8 <sup>th</sup> , 2014	<b>CHANGES TO USER GUIDE</b> - The kit content was not changed.	
	The extra chloroform extraction was removed. Workflow, preparation table, number of PLG-tubes, volumes of user-supplied reagents and the phenol-chloroform extraction protocol were adapted accordingly.	5, 6, 7, 12
	Storage of all kit components can now be at +2 to +8 °C (+4 °C).	6
	Incubation and centrifugation times were shortened.	12 - 17
	Isopropanol volume increased to 1.75 x to maximize miRNA recovery.	13
	Max. loading volume of purification column increased from 600 µl to 800 µl.	13, 15, 16
	No re-elution but optional second elution into new micro-tube.	14, 17
	Short Procedures were adapted accordingly.	19 - 20
August 27 <sup>th</sup> , 2014	<b>CHANGES TO USER GUIDE</b>	
	Figure 3 updated to include miRNA spike-in experiment.	21
September 9 <sup>th</sup> , 2014	<b>FORMAT CHANGES TO REFERENCE CARD</b>	
May 8 <sup>th</sup> , 2015	<b>CHANGES TO USER GUIDE</b>	
	General text changes to account for consistency.	
	Note added on RNasin 230 nm absorption to be considered for OD blanking.	14, 17
	Optional RNasin addition added to Short Procedures.	18, 19

A decorative graphic consisting of several translucent blue spheres of various sizes connected by thin, light blue lines, creating a network-like structure across the page.

## SPLIT RNA Extraction Kit · User Guide

Lexogen GmbH  
Campus Vienna Biocenter 5  
1030 Vienna, Austria  
Telephone: +43 (0) 1 345 1212-41  
Fax: +43 (0) 1 345 1212-99  
E-mail: [info@lexogen.com](mailto:info@lexogen.com)  
© Lexogen GmbH, 2015