



# User Information - DNA Removal Add-on

Residual DNA should be removed from FFPE RNA prior to rRNA depletion. Most FFPE RNA extraction kits contain a DNase I treatment, however, this may not be complete, depending on the amount of DNA that was co-purified. Lexogen offers a DNA Removal Add-on that can be seamlessly integrated upstream of RiboCop HMR V2 rRNA depletion. There is no need to purify samples between DNA removal and rRNA depletion. Total FFPE RNA from 5 ng - 1 µg can be used as input for DNA removal.

**ATTENTION:** Before starting this protocol, please read the [General Guidelines for Lexogen Kits](#), which are available online. These provide a detailed overview of RNA and kit component handling, as well as general RNA input requirements.

## 1. Kit Components

### DNA Removal Add-on

Kit Component	Tube Label	Volume*			Storage
		24 preps	96 preps	384 preps	
10x DNA Removal Buffer	DRB ●	27 µl	106 µl	4x 106 µl	 -20 °C
DNA Removal Enzyme	DRE ●	27 µl	106 µl	4x 106 µl	 -20 °C

\*including ≥10 % surplus

## 2. User-Supplied Consumables and Equipment

### Reagents

- RNase-free water.

### Equipment

- Benchtop centrifuge (12,000 x g, rotor compatible with 1.5 ml tubes or 3,000 x g, rotor compatible with 96-well plates).
- Calibrated single-channel and multi-channel pipettes for handling 1 µl to 1,000 µl volumes.
- Thermocycler.
- UV-spectrophotometer to quantify RNA.
- Ice bath, ice box, ice pellets, or benchtop cooler (-20 °C for enzymes).

## Labware

- Suitable low-binding pipette tips (pipette tips with aerosol barriers recommended).
- 1.5 ml reaction tubes, low binding, certified ribonuclease-free.
- 200 µl PCR tubes or 96-well plates and caps or sealing foil.
- Vortex mixer.

## Optional Equipment

- Automated microfluidic electrophoresis station (e.g., Agilent Technologies 2100 Bioanalyzer).
- Benchtop fluorometer and appropriate assays (for RNA quality control).
- Agarose gels, dyes, and electrophoresis rig (for RNA quality control).

## 3. Short Protocol

DNA Removal	
<b>DRB ●</b>	– thawed at RT
<b>DRE ●</b>	– keep on ice or at -20 °C
Thermocycler	37 °C, 10 min

Up to 8 µl total FFPE RNA, ranging from 5 ng - 1 µg can be inserted into the DNA removal step.

**ATTENTION:** For the preparation of mastermixes include a 10 % surplus per reaction.

**EXAMPLE:** Step 1 for 24 preps: 26.4 µl **DRB ●** (= 1 µl x 24 rxn x 1.1)  
+ 26.4 µl **DRE ●** (= 1 µl x 24 rxn x 1.1)

resulting in a total of 52.8 µl, which is sufficient for multi-channel pipetting.

1 Prepare a mastermix of 1 µl 10x DNA Removal Buffer (**DRB ●**) and 1 µl DNA Removal Enzyme (**DRE ●**) per sample. Mix thoroughly and spin down briefly.

2 Add 2 µl of the **DRB / DRE** mastermix to 8 µl of FFPE RNA sample. If a smaller volume of RNA is used, add RNase-free water to a total reaction volume of 10 µl. Mix thoroughly and quickly spin down.

3 Incubate for 10 minutes at 37 °C.  Safe stopping point. After completing the reaction, the samples can be stored at -20 °C.

Add 16 µl RNase-free water to a total volume of 26 µl and insert the sample into step 2 of RiboCop rRNA depletion and continue with the protocol for CORALL FFPE RNA-Seq as described in 219UG781.

4 Alternatively, purify the reaction using spin-column purification or magnetic beads suitable for RNA purification before proceeding to downstream reactions. For more information, please contact [support@lexogen.com](mailto:support@lexogen.com).